

Computerized pressure AlgoMeter for assessing deep tissue pain perception in Muscle Pain Syndromes & Fibromyalgia

Software-based REAL TIME visual & auditory feedback to control & monitor Applied Pressure Rates

- Real-time visual & auditory feedback
- Available as Stand-alone product or add-on to PATHWAY System
- Comprehensive test statistics
- Body site selection
- USB Patient Response Unit
- Customizable outcome-based color graphic reports
- Calibration verification weight
- Ramp & Hold Protocol
- Normative data presentation



Hardware Advantages

- Ergonomic handle for easy grip
- Rechargeable 50 hour battery
- Support for various tip sizes




Patient Response Unit

Test Management

- Graphic display of test including applied pressure change rate
- Display test statistics according to selected test methods

Software Capabilities

- Allows managing patients, test paradigms and results
- Pressure units support includes KPa, Kg/cm² and lbs/in²
- Patient response to pressure threshold and tolerance tests can be recorded
- Body site selection
- Test results can be mapped according to a predefined body site
- Test Results can be saved, exported to Excel and printed as a customizable report



AlgoMeter test screen

Medoc Ltd. | Headquarters

1 Ha'Dekel St. Ramat Yishai, Israel 30095
phone: +972-4-903.8800 | info@medoc-web.com

Medoc USA | Compass Medical Technologies, Inc

1502 West Highway 54, Durham, North Carolina 27707
phone: 919.402.9600 | medoc@mindspring.com

www.medoc-web.com